

Product Information

Product SKU:	CAB12883	Gene ID:	56985	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Mouse,Rat

Additional Information

Observed MW:	39kDa	Conjugate:	Unconjugated
Calculated MW:	40kDa	Isotype:	IgG

Immunogen Information

**Background:** Predicted to enable 2',3'-cyclic-nucleotide 2'-phosphodiesterase activity; manganese ion binding activity; and pyrophosphatase activity. Predicted to be located in cytosol.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	MDS006; C17orf48; NBLA03831; ADPRM
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 203-342 of human ADPRM (NP_064618.3).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.