

# PARP8 Rabbit Polyclonal Antibody



CAB12894

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

96kDa

### Calculated MW:

91kDa/95kDa

### Applications:

WB IHC

### Reactivity:

Human, Mouse, Rat

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000 IHC 1:50  
- 1:200

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

### Immunogen information

#### Gene ID:

79668

#### Uniprot

Q8N3A8

#### Synonyms:

PARP8; ARTD16; pART16

#### Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human PARP8 (NP\_001171527.1).

#### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---

