

Product Information

Product SKU:	CAB12911	Gene ID:	84303	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	26kDa	Conjugate:	Unconjugated
Calculated MW:	26kDa	Isotype:	IgG

Immunogen Information

Background: Involved in cellular response to DNA damage stimulus and cristae formation. Located in cytosol and mitochondrial inner membrane. Part of MICOS complex.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	CHCM1; Mic25; MICOS25; PPP1R23; CHCHD6
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-235 of human CHCHD6 (NP_115719.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.