

SFXN1 Rabbit Polyclonal Antibody

CAB12954



Product Information

Product SKU:	CAB12954	Gene ID:	94081	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	36kDa	Conjugate:	Unconjugated
Calculated MW:	36kDa	Isotype:	IgG

Immunogen Information

Background:	Enables D-serine transmembrane transporter activity and L-serine transmembrane transporter activity. Involved in D-serine transport; L-serine transport; and serine import into mitochondrion. Located in mitochondrion. Is integral component of mitochondrial inner membrane.
Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	TCC; SLC56A1; SFXN1
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 87-322 of human SFXN1 (NP_073591.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.