

WSB2 Rabbit Polyclonal Antibody



CAB12970

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

40-43kDa

Calculated MW:

21kDa/45kDa/47kDa

Applications:

WB IHC

Reactivity:

Mouse, Rat

Protein Background

This gene encodes a member of the WD-protein subfamily. The encoded protein contains five WD-repeats spanning most of the protein and an SOCS box in the C-terminus. The SOCS box may act as a bridge between specific substrate-binding domains and E3 ubiquitin protein ligases. Three transcript variants encoding different isoforms have been found for this gene.

Immunogen information

Gene ID:

55884

Uniprot

Q9NYS7

Synonyms:

WSB2; SBA2

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000 IHC 1:50
- 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human WSB2 (NP_061109.1).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images

