
Product Information

| | | | | | |
|---------------------|----------|----------------------|--------|--------------------|-------------|
| Product SKU: | CAB12970 | Gene ID: | 55884 | Size: | 20uL, 100uL |
| Clone No: | - | Host Species: | Rabbit | Reactivity: | Mouse,Rat |

Additional Information

| | | | |
|-----------------------|----------|-------------------|--------------|
| Observed MW: | 40-43kDa | Conjugate: | Unconjugated |
| Calculated MW: | 45kDa | Isotype: | IgG |

Immunogen Information

Background: This gene encodes a member of the WD-protein subfamily. The encoded protein contains five WD-repeats spanning most of the protein and an SOCS box in the C-terminus. The SOCS box may act as a bridge between specific substrate-binding domains and E3 ubiquitin protein ligases. Three transcript variants encoding different isoforms have been found for this gene.

Recommended Dilution: WB,1:500 - 1:2000 IHC-P,1:50 - 1:200

Synonyms: SBA2; WSB2

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human WSB2 (NP_061109.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.