

Product Information

Product SKU:	CAB13006	Gene ID:	51063	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Rat

Additional Information

Observed MW:	36kDa	Conjugate:	Unconjugated
Calculated MW:	36kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable cation channel activity. Involved in positive regulation of apoptotic process.
Predicted to be integral component of plasma membrane.

Recommended Dilution:	WB,1:100 - 1:500
Synonyms:	FAM26B; CALHM2
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-323 of human CALHM2 (NP_057000.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.