

CAB13073

Product Information

Product SKU:	CAB13073	Gene ID:	10906	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	78kDa	Conjugate:	Unconjugated
Calculated MW:	65kDa	Isotype:	IgG

Immunogen Information

Background: The innate immune system confers host defense against viral and microbial infection, and TRAFD1 is a negative feedback regulator that controls excessive immune responses (Sanada et al., 2008 [PubMed 18849341]).

Recommended Dilution:	WB,1:500 - 1:1000
Synonyms:	FLN29; TRAFD1
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human TRAFD1 (NP_006691.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.