
Product Information

Product SKU:	CAB13182	Gene ID:	6999	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	48kDa	Conjugate:	Unconjugated
Calculated MW:	48kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes a heme enzyme that plays a critical role in tryptophan metabolism by catalyzing the first and rate-limiting step of the kynurenine pathway. Increased activity of the encoded protein and subsequent kynurenine production may also play a role in cancer through the suppression of antitumor immune responses, and single nucleotide polymorphisms in this gene may be associated with autism.

Recommended Dilution: WB,1:500 - 1:1000 IF/ICC,1:100 - 1:500

Synonyms: TO; TDO; TPH2; TRPO; HYPTRP; TDO2

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 100-200 of human TDO2 (NP_005642.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.