
Product Information

Product SKU:	CAB13185	Gene ID:	7763	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Mouse

Additional Information

Observed MW:	23kDa	Conjugate:	Unconjugated
Calculated MW:	23kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable DNA binding activity and zinc ion binding activity. Predicted to act upstream of or within several processes, including face development; fibroblast migration; and platelet-derived growth factor receptor signaling pathway. Predicted to be located in cytoplasm.

Recommended Dilution: WB, 1:500 - 1:1000

Synonyms: ZA20D2; ZNF216; ZFAND5A; ZFAND5

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 40-150 of human ZFAND5 (NP_001095890.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.