

Product Information

Product SKU:	CAB13191	Gene ID:	29889	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	84kDa	Conjugate:	Unconjugated
Calculated MW:	84kDa	Isotype:	IgG

Immunogen Information

Background: Enables RNA binding activity. Predicted to be involved in ribosome biogenesis. Located in nucleolus.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	NGP1; Nog2; Nug2; Ngp-1; HUMAUANTIG; GNL2
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human GNL2 (NP_037417.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.