

OAF Rabbit Polyclonal Antibody

CAB13202



Product Information

Product SKU:	CAB13202	Gene ID:	220323	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Mouse,Rat

Additional Information

Observed MW:	30kDa	Conjugate:	Unconjugated
Calculated MW:	31kDa	Isotype:	IgG

Immunogen Information

**Background:** OAF (Out At First Homolog) is a Protein Coding gene. Diseases associated with OAF include Spondylocarpotarsal Synostosis Syndrome and Excessive Tearing.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	NS5ATP13TP2; OAF
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 28-150 of human OAF (NP_848602.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.