

# TESC Rabbit Polyclonal Antibody



CAB13219

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

Refer to Figures

### Calculated MW:

18kDa/21kDa/24kDa

### Applications:

WB

### Reactivity:

Mouse

## Antibody Information

### Recommended dilutions:

WB 1:1000 - 1:2000

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

### Immunogen information

#### Gene ID:

54997

#### Uniprot

Q96BS2

#### Synonyms:

TESC; CHP3; TSC; tescalcin

#### Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 135-214 of human TESC (NP\_060369.3).

#### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---