

Product Information

Product SKU:	CAB13234	Gene ID:	115761	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Rat

Additional Information

Observed MW:	21kDa	Conjugate:	Unconjugated
Calculated MW:	21kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes a tumor suppressor related to the ADP-ribosylation factor (ARF) family of proteins. The encoded protein may play a role in apoptosis in a caspase-dependent manner. Polymorphisms in this gene have been associated with some familial cancers.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	ARLTS1; ARL11
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-196 of human ARL11 (NP_612459.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.