

TBCB Rabbit Polyclonal Antibody



CAB13248

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

37kDa

Calculated MW:

21kDa/27kDa

Applications:

WB IHC IF

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:1000 - 1:2000 IHC
1:50 - 1:200 IF 1:50 - 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

1155

Uniprot

Q99426

Synonyms:

TBCB; CG22; CKAP1; CKAPI



Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 1-244 of human TBCB (NP_001272.2).


Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.


Product Images



Immunohistochemistry of paraffin-embedded rat brain using TBCB antibody (CAB13248) at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded human placenta using TBCB antibody (CAB13248) at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded mouse kidney using TBCB antibody (CAB13248) at dilution of 1:100 (40x lens).