

# CD164 Rabbit Polyclonal Antibody



CAB13254

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

Refer to Figures

### Calculated MW:

16kDa/19kDa/20kDa

### Applications:

WB

### Reactivity:

Human

## Antibody Information

### Recommended dilutions:

WB 1:1000 - 1:2000

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

### Immunogen information

#### Gene ID:

8763

#### Uniprot

Q04900

#### Synonyms:

CD164; DFNA66; MGC-24; MGC-24v; MUC-24; endolyn

#### Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 24-162 of human CD164 (NP\_006007.2).

#### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---