

# USP6 Rabbit Polyclonal Antibody



CAB13255

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

Refer to Figures

### Calculated MW:

89kDa/121kDa/158kDa

### Applications:

WB

### Reactivity:

Human

## Antibody Information

### Recommended dilutions:

WB 1:1000 - 1:2000

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

### Immunogen information

#### Gene ID:

9098

#### Uniprot

P35125

#### Synonyms:

USP6; HRP1; TRE17; TRE2; TRESMCR; Tre-2; USP6-short

#### Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 55-120 of human USP6 (NP\_004496.2).

#### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---