

YTHDF1 Rabbit Polyclonal Antibody



CAB13260

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

61kDa

Calculated MW:

20kDa/60kDa

Applications:

WB IF

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000 IF 1:50 - 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

54915

Uniprot

Q9BYJ9

Synonyms:

YTHDF1; C20orf21


Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human YTHDF1 (NP_060268.2).


Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.


Product Images



Western blot analysis of extracts of various cell lines, using YTHDF1 antibody (CAB13260) at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (CABS014) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit (CABM00020). Exposure time: 60s.



Immunofluorescence analysis of L929 cells using YTHDF1 Rabbit pAb (CAB13260) at dilution of 1:100. Blue: DAPI for nuclear staining.



Immunofluorescence analysis of U2OS cells using YTHDF1 Rabbit pAb (CAB13260) at dilution of 1:100. Blue: DAPI for nuclear staining.