

VGLL3 Rabbit Polyclonal Antibody



CAB13263

Product Information

Size:

50uL, 100uL, 200uL

Observed MW:

Calculated MW:

35kDa/36kDa

Applications:

IHC

Reactivity:

Human, Mouse

Antibody Information

Recommended dilutions:

IHC 1:50 - 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

389136

Uniprot

A8MV65

Synonyms:

VGLL3; VGL-3; VGL3

Immunogen:

A synthetic peptide of human VGLL3

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
