
Product Information

Product SKU:	CAB13282	Gene ID:	716	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	88kDa	Conjugate:	Unconjugated
Calculated MW:	77kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes a serine protease, which is a major constituent of the human complement subcomponent C1. C1s associates with two other complement components C1r and C1q in order to yield the first component of the serum complement system. Defects in this gene are the cause of selective C1s deficiency.

Recommended Dilution: WB, 1:500 - 1:2000

Synonyms: EDSPD2; C1S

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 150-250 of human C1S (NP_001725.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.