
Product Information

Product SKU:	CAB13411	Gene ID:	26019	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Rat

Additional Information

Observed MW:	170kDa	Conjugate:	Unconjugated
Calculated MW:	148kDa	Isotype:	IgG

Immunogen Information

Background:	This gene encodes a protein that is part of a post-splicing multiprotein complex involved in both mRNA nuclear export and mRNA surveillance. mRNA surveillance detects exported mRNAs with truncated open reading frames and initiates nonsense-mediated mRNA decay (NMD). When translation ends upstream from the last exon-exon junction, this triggers NMD to degrade mRNAs containing premature stop codons. This protein is located in the perinuclear area. It interacts with translation release factors and the proteins that are functional homologs of yeast Upf1p and Upf3p. Two splice variants have been found for this gene; both variants encode the same protein.
Recommended Dilution:	WB,1:500 - 1:2000 IHC-P,1:50 - 1:200 IP,0.5µg-4µg antibody for 200µg-400µg extracts of whole cells
Synonyms:	HUPF2; RENT2; smg-3; UPF2
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1153-1272 of human UPF2 (NP_542166.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.