

NAT8B Rabbit Polyclonal Antibody



CAB13428

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Calculated MW:

25kDa

Applications:

IF

Reactivity:

Human

Antibody Information

Recommended dilutions:

IF 1:50 - 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

51471

Uniprot

Q9UHF3

Synonyms:

NAT8B; CML2; Hcml2; NAT8BP


Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 75-165 of human NAT8B (NP_057431.2).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunofluorescence analysis of A549 cells using NAT8B antibody (CAB13428). Blue: DAPI for nuclear staining.