

# TIMM17A Rabbit Polyclonal Antibody



CAB13578

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

18kDa

### Calculated MW:

18kDa

### Applications:

WB IHC IF

### Reactivity:

Human, Mouse

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000 IHC 1:50  
- 1:200 IF 1:50 - 1:200

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

### Immunogen information

#### Gene ID:

10440

#### Uniprot

Q99595

#### Synonyms:

TIMM17A; TIM17; TIM17A

#### Immunogen:


Recombinant fusion protein containing a sequence corresponding to amino acids 1-171 of human TIMM17A (NP\_006326.1).

#### Storage:


Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images


---



Western blot analysis of extracts of various cell lines, using TIMM17A antibody (CAB13578) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (CABS014) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit (CABM00020). Exposure time: 30s.



Immunohistochemistry of paraffin-embedded human liver cancer using TIMM17A antibody (CAB13578) at dilution of 1:100 (40x lens).



Immunofluorescence analysis of U2OS cells using TIMM17A antibody (CAB13578). Blue: DAPI for nuclear staining.