

CAB13582

Product Information

Product SKU:	CAB13582	Gene ID:	11091	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	37kDa	Conjugate:	Unconjugated
Calculated MW:	37kDa	Isotype:	IgG

Immunogen Information

Background:	This gene encodes a member of the WD repeat protein family. WD repeats are minimally conserved regions of approximately 40 amino acids typically bracketed by gly-his and trp-aspartate (GH-WD), which may facilitate formation of heterotrimeric or multiprotein complexes. Members of this family are involved in a variety of cellular processes, including cell cycle progression, signal transduction, apoptosis, and gene regulation. This protein contains 7 WD repeats. Alternatively spliced transcript variants encoding the same protein have been identified.
Recommended Dilution:	WB, 1:500 - 1:2000
Synonyms:	BIG3; SWD3; BIG-3; CFAP89; WDR5
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-334 of human WDR5 (NP_438172.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.