

RAD54L2 Rabbit Polyclonal Antibody



CAB13585

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Calculated MW:

162kDa

Applications:

IF

Reactivity:

Human

Antibody Information

Recommended dilutions:

IF 1:50 - 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

23132

Uniprot

Q9Y4B4

Synonyms:

RAD54L2; ARIP4; HSPC325; SRISNF2L


Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 1178-1467 of human RAD54L2 (NP_055921.2).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunofluorescence analysis of MCF-7 cells using RAD54L2 antibody (CAB13585). Blue: DAPI for nuclear staining.