
Product Information

Product SKU:	CAB13585	Gene ID:	23132	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Rat

Additional Information

Observed MW:	-	Conjugate:	Unconjugated
Calculated MW:	163kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable ATP hydrolysis activity; protein kinase binding activity; and transcription coregulator activity. Predicted to be involved in DNA duplex unwinding. Predicted to act upstream of or within positive regulation of transcription by RNA polymerase II. Predicted to be located in nuclear speck. Predicted to be active in nucleus.

Recommended Dilution: IHC-P,1:50 - 1:200 IF/ICC,1:50 - 1:200

Synonyms: ARIP4; HSPC325; SRISNF2L; RAD54L2

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1178-1467 of human RAD54L2 (NP_055921.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.