

P4HTM Rabbit Polyclonal Antibody



CAB13593

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Calculated MW:

47kDa/56kDa/63kDa

Applications:

IF

Reactivity:

Human

Antibody Information

Recommended dilutions:

IF 1:50 - 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

The product of this gene belongs to the family of prolyl 4-hydroxylases. This protein is a prolyl hydroxylase that may be involved in the degradation of hypoxia-inducible transcription factors under normoxia. It plays a role in adaptation to hypoxia and may be related to cellular oxygen sensing. Alternatively spliced variants encoding different isoforms have been identified.

Immunogen information

Gene ID:

54681

Uniprot

Q9NXG6

Synonyms:

P4HTM; EGLN4; HIFPH4; P4H-TM; PH-4; PH4; PHD4


Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 294-563 of human P4HTM (NP_808807.2).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunofluorescence analysis of HeLa cells using P4HTM antibody (CAB13593). Blue: DAPI for nuclear staining.