

# CCAR1 Rabbit Polyclonal Antibody



CAB13595

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

132KDa

### Calculated MW:

131kDa/132kDa

### Applications:

WB IF IP

### Reactivity:

Human, Mouse, Rat

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000 IF 1:50 -  
1:200 IP 1:50 - 1:100

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

### Immunogen information

#### Gene ID:

55749

#### Uniprot

Q8IX12

#### Synonyms:

CCAR1

#### Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human CCAR1 (NP\_060707.2).

#### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

