

Product Information

| | | | | | |
|--------------|----------|---------------|--------|-------------|-------------|
| Product SKU: | CAB13649 | Gene ID: | 5934 | Size: | 20uL, 100uL |
| Clone No: | - | Host Species: | Rabbit | Reactivity: | Human,Mouse |

Additional Information

| | | | |
|----------------|--------|------------|--------------|
| Observed MW: | 130kDa | Conjugate: | Unconjugated |
| Calculated MW: | 128kDa | Isotype: | IgG |

Immunogen Information

Background: Enables promoter-specific chromatin binding activity. Involved in regulation of lipid kinase activity. Acts upstream of or within negative regulation of gene expression. Located in chromosome; cytosol; and nuclear lumen.

| | |
|-----------------------|---|
| Recommended Dilution: | WB,1:500 - 1:1000 IF/ICC,1:50 - 1:200 |
| Synonyms: | Rb2; P130; BRUWAG; RBL2 |
| Purification Method: | Affinity purification |
| Immunogen: | A synthetic peptide corresponding to a sequence within amino acids 413-613 of human RBL2 (NP_005602.3). |
| Storage: | Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3. |