

# RBL2 Rabbit Polyclonal Antibody



CAB13649

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

Refer to Figures

### Calculated MW:

59kDa/128kDa

### Applications:

WB

### Reactivity:

Human

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

### Immunogen information

#### Gene ID:

5934

#### Uniprot

Q08999

#### Synonyms:

RBL2; P130; Rb2

#### Immunogen:

A synthetic peptide corresponding to a sequence within amino acids 100-200 of human RBL2 (NP\_005602.3).

#### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---