

KIAA1456 Rabbit Polyclonal Antibody



CAB13657

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Calculated MW:

Applications:

IF

Reactivity:

Human

Antibody Information

Recommended dilutions:

IF 1:50 - 1:100

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

57604

Uniprot

Q9P272

Synonyms:

KIAA1456; C8orf79; TRM9L

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 110-330 of human KIAA1456 (NP_065895.2).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
