

Product Information

Product SKU:	CAB13692	Gene ID:	375189	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Mouse,Rat

Additional Information

Observed MW:	14kDa	Conjugate:	Unconjugated
Calculated MW:	14kDa	Isotype:	IgG

Immunogen Information

**Background:** Predicted to enable actin monomer binding activity. Predicted to be involved in sequestering of actin monomers. Predicted to be located in cytoskeleton. Predicted to be active in cell cortex.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	PFN4
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-129 of human PFN4 (NP_955378.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.