

Product Information

Product SKU:	CAB13716	Gene ID:	27071	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	28kDa	Conjugate:	Unconjugated
Calculated MW:	32kDa	Isotype:	IgG

Immunogen Information

Background: Enables phosphatidylinositol-3,4,5-trisphosphate binding activity and phosphatidylinositol-3,4-bisphosphate binding activity. Predicted to be involved in signal transduction. Located in plasma membrane.

Recommended Dilution:	WB,1:500 - 1:2000 IF/ICC,1:50 - 1:200
Synonyms:	BAM32; DAPP1
Purification Method:	Affinity purification
Immunogen:	A synthetic peptide corresponding to a sequence within amino acids 100-273 of human DAPP1 (NP_055210.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.