

ASCC2 Rabbit Polyclonal Antibody

CAB13789



Product Information

Product SKU:	CAB13789	Gene ID:	84164	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	110kDa	Conjugate:	Unconjugated
Calculated MW:	86kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable ubiquitin binding activity. Involved in regulation of transcription, DNA-templated; rescue of stalled ribosome; and ribosome-associated ubiquitin-dependent protein catabolic process. Located in nucleus. Part of activating signal cointegrator 1 complex.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	p100; ASC1p100; ASCC2
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 488-757 of human ASCC2 (NP_115580.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.