

Product Information

Product SKU:	CAB13791	Gene ID:	9610	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Mouse,Rat

Additional Information

Observed MW:	84kDa	Conjugate:	Unconjugated
Calculated MW:	84kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable small GTPase binding activity. Predicted to be involved in endocytosis; regulation of catalytic activity; and signal transduction. Predicted to act upstream of or within associative learning; memory; and negative regulation of synaptic plasticity. Located in cytoplasm and plasma membrane.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	RIN1
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human RIN1 (NP_004283.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.