

RBSN Rabbit Polyclonal Antibody



CAB13792

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

23kDa/88kDa

Applications:

WB

Reactivity:

Mouse

Protein Background

This gene encodes a protein that belongs to the FYVE zinc finger family of proteins. The encoded protein interacts with Ras-related proteins that regulate membrane trafficking. A missense mutation in this gene is associated with a defect in the early endocytic pathway. Alternative splicing results in multiple transcript variants.

Immunogen information

Gene ID:

64145

Uniprot

Q9H1K0

Synonyms:

RBSN; Rabenosyn-5; ZFYVE20

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human RBSN (NP_071735.2).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
