

ZNF227 Rabbit Polyclonal Antibody



CAB13793

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

85kDa/92kDa

Applications:

WB IHC

Reactivity:

Human

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000 IHC 1:50
- 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

7770

Uniprot

Q86WZ6

Synonyms:

ZNF227

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 60-230 of human ZNF227 (NP_872296.1).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
