

CAB13847

Product Information

Product SKU:	CAB13847	Gene ID:	79866	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	90kDa	Conjugate:	Unconjugated
Calculated MW:	61kDa	Isotype:	IgG

Immunogen Information

**Background:** BORA is an activator of the protein kinase Aurora A (AURKA; MIM 603072), which is required for centrosome maturation, spindle assembly, and asymmetric protein localization during mitosis (Hutterer et al., 2006 [PubMed 16890155]).

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	C13orf34; BORA
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 320-619 of human BORA (NP_079084.3).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.