

IGFBP5 Rabbit Polyclonal Antibody



CAB13858

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Calculated MW:

30kDa

Applications:

WB IF

Reactivity:

Human, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000 IF 1:50 - 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

3488

Uniprot

P24593

Synonyms:

IGFBP5; IBP5


Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 21-272 of human IGFBP5 (NP_000590.1).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunofluorescence analysis of HeLa cells using IGFBP5 antibody (CAB13858).