

# MGAT5B Rabbit Polyclonal Antibody



CAB13870

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

89kDa

### Calculated MW:

89kDa/90kDa

### Applications:

WB

### Reactivity:

Mouse

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

### Immunogen information

#### Gene ID:

146664

#### Uniprot

Q3V5L5

#### Synonyms:

MGAT5B; GnT-IX; GnT-VB; alpha-1

#### Immunogen:


Recombinant fusion protein containing a sequence corresponding to amino acids 616-790 of human MGAT5B (NP\_653278.2).

#### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---



Western blot analysis of extracts of mouse brain, using MGAT5B antibody (CAB13870) at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (CABS014) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit (CABM00021). Exposure time: 90s.