

SPSB1 Rabbit Polyclonal Antibody

CAB13886



Product Information

| | | | | | |
|---------------------|----------|----------------------|--------|--------------------|-----------------|
| Product SKU: | CAB13886 | Gene ID: | 80176 | Size: | 20uL, 100uL |
| Clone No: | - | Host Species: | Rabbit | Reactivity: | Human,Mouse,Rat |

Additional Information

| | | | |
|-----------------------|-------|-------------------|--------------|
| Observed MW: | 31kDa | Conjugate: | Unconjugated |
| Calculated MW: | 31kDa | Isotype: | IgG |

Immunogen Information

Background: Enables ubiquitin ligase-substrate adaptor activity. Involved in protein ubiquitination and ubiquitin-dependent protein catabolic process. Located in cytosol.

| | |
|------------------------------|--|
| Recommended Dilution: | WB,1:500 - 1:2000 |
| Synonyms: | SSB1; SSB-1; SPSB1 |
| Purification Method: | Affinity purification |
| Immunogen: | Recombinant fusion protein containing a sequence corresponding to amino acids 1-60 of human SPSB1 (NP_079382.2). |
| Storage: | Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3. |