

Product Information

Product SKU:	CAB13888	Gene ID:	140735	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	12kDa	Conjugate:	Unconjugated
Calculated MW:	10kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable dynein intermediate chain binding activity and dynein light intermediate chain binding activity. Predicted to be involved in cilium assembly. Located in 9+0 non-motile cilium and centrosome. Is active in glutamatergic synapse and postsynapse.

Recommended Dilution:	WB,1:500 - 1:2000 IF/ICC,1:50 - 1:200
Synonyms:	Dlc2; DNCL1B; RSPH22; DYNLL2
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-89 of human DYNLL2 (NP_542408.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.