

CAB13893

Product Information

Product SKU:	CAB13893	Gene ID:	83714	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	31kDa	Conjugate:	Unconjugated
Calculated MW:	31kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable aspartic-type endopeptidase activity. Predicted to be involved in proteolysis. Predicted to act upstream of or within Notch signaling pathway and negative regulation of transcription by RNA polymerase II. Located in cytoplasm.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	NRIP2
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 140-281 of human NRIP2 (NP_113662.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.