

SLC5A6 Rabbit Polyclonal Antibody



CAB13904

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Calculated MW:

68kDa

Applications:

IF

Reactivity:

Human

Antibody Information

Recommended dilutions:

IF 1:50 - 1:100

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

8884

Uniprot

Q9Y289

Synonyms:

SLC5A6; SMVT

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 546-635 of human SLC5A6 (NP_066918.2).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
