

CAB13914

Product Information

Product SKU:	CAB13914	Gene ID:	79004	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Rat

Additional Information

Observed MW:	32kDa	Conjugate:	Unconjugated
Calculated MW:	32kDa	Isotype:	IgG

Immunogen Information

**Background:** Predicted to enable ubiquitin binding activity. Acts upstream of or within negative regulation of cytokine production involved in inflammatory response and negative regulation of macrophage cytokine production. Located in cytosol; nuclear membrane; and nucleoplasm.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	C10orf66; bA18I14.5; CUEDC2
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-150 of human CUEDC2 (NP_076945.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.