

Product Information

Product SKU:	CAB1404	Gene ID:	1903	Size:	20uL, 100uL
Clone No:	ARC1877	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	42kDa	Conjugate:	Unconjugated
Calculated MW:	42kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes a member of the EDG family of receptors, which are G protein-coupled receptors. This protein has been identified as a functional receptor for sphingosine 1-phosphate and likely contributes to the regulation of angiogenesis and vascular endothelial cell function.

Recommended Dilution:	WB,1:500 - 1:1000
Synonyms:	EDG3; LPB3; S1P3; EDG-3; C9orf47; C9orf108; bA791O21.3; S1PR3
Purification Method:	Affinity purification
Immunogen:	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human S1PR3 (Q99500).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.