

RAD51B Rabbit Polyclonal Antibody



CAB14056

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Calculated MW:

25kDa/38kDa/40kDa/42kDa

Applications:

IF

Reactivity:

Human

Antibody Information

Recommended dilutions:

IF 1:50 - 1:100

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

The protein encoded by this gene is a member of the RAD51 protein family. RAD51 family members are evolutionarily conserved proteins essential for DNA repair by homologous recombination. This protein has been shown to form a stable heterodimer with the family member RAD51C, which further interacts with the other family members, such as RAD51, XRCC2, and XRCC3. Overexpression of this gene was found to cause cell cycle G1 delay and cell apoptosis, which suggested a role of this protein in sensing DNA damage. Rearrangements between this locus and high mobility group AT-hook 2 (HMGA2, GenelD 8091) have been observed in uterine leiomyomata.

Immunogen information

Gene ID:

5890

Uniprot

O15315

Synonyms:

RAD51B; R51H2; RAD51L1; REC2

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 125-384 of human RAD51B (NP_598193.2).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
