

N4BP1 Rabbit Polyclonal Antibody



CAB14094

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

120kDa

Calculated MW:

99kDa

Applications:

WB

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

80750

Uniprot

Synonyms:

N4BP1


Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 10-343 of mouse N4BP1 (NP_085040.2).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot - N4BP1 Rabbit pAb (CAB14094)