

# FAIM2 Rabbit Polyclonal Antibody



CAB14110

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

Refer to figures

### Calculated MW:

30kDa/35kDa

### Applications:

WB

### Reactivity:

Mouse

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

### Immunogen information

#### Gene ID:

23017

#### Uniprot

Q9BWQ8

#### Synonyms:

FAIM2; LFG; LFG2; NGP35; NMP35; TMBIM2

#### Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 1-105 of human FAIM2 (NP\_036438.2).

#### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---