

G0S2 Rabbit Polyclonal Antibody

CAB14125



Product Information

Product SKU:	CAB14125	Gene ID:	50486	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	11kDa	Conjugate:	Unconjugated
Calculated MW:	11kDa	Isotype:	IgG

Immunogen Information

Background: Involved in extrinsic apoptotic signaling pathway and positive regulation of extrinsic apoptotic signaling pathway. Located in mitochondrion.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	G0S2
Purification Method:	Affinity purification
Immunogen:	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human G0S2 (NP_056529.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.